


Trends on Montel X-ray Optics for Inelastic Scattering and Pinholes for Synchrotron Beamlines

F. Hertlein, U. Heidorn, C. Umland, C. Michaelsen, J. Wiesmann

Incoatec GmbH, Max-Planck-Str. 2, 21502 Geesthacht, Germany

Montel Optics for Synchrotron Applications in Different Sizes

Montel optics consist of two mirrors mounted side-by-side in an L-shape arrangement enabling a 2-dimensional beam shaping. A Montel optics with two elliptically shaped mirrors is point focusing, whereas two parabolic mirrors enable a collimated beam. A line focus is created with a hybrid optics, a combination of an elliptic and a parabolic mirror. High quality multilayer optics are essential for an excellent beam shaping with homogeneous beam properties. The Montel optics accumulate a lot of flux within a well-shaped, gaussian-like spot of expected size measured by 2D detectors or pin diodes. Nowadays, Montel optics are also used at synchrotrons, where they substitute the KB (Kirkpatrick-Baez) mirrors enabling a more compact design.




Montel Optics - 100-250 mm in length
Different cross sections from 40 x 40 mm to 10 x 10 mm. First optics, with slope errors < 2 arcsec, were sold to NSLS, DLS, APS and ESRF. They are used at inelastic scattering beamlines.

Multi-stripe Multilayer Optics

At imaging beamlines multilayer optics are often used as double crystal multilayer monochromators (DCMM). For example, tomography needs a homogeneous and stable beam profile, in order to perform optimal background corrections. Because of the high coherence of the radiation, the optical components must be designed with particular care in order to avoid a deterioration of the beam quality. Multilayer coatings with up to 5 stripes were produced with films homogeneities < 0.2% as well as with lateral gradients.

Stripe A: [Ru/C]100, $d=40 \text{ \AA}$, $\gamma=0.5$,
 $R > 80\%$ for $10 < E < 22 \text{ keV}$
Midspace: Si<111>, roughness 0.1 nm,
slope error 0.04"
Stripe B: [W/Si]100, $d=30 \text{ \AA}$, $\gamma=0.5$,
 $R > 80\%$ for $22 < E < 45 \text{ keV}$



Three-striped multilayer optics for tomographic microscopy and coherent radiology, with an optimized coating for different beam energies (TOMCAT at SLS, Switzerland, Data courtesy of M. Stampantonii).

Si, orientation <0.05° 300 x 55 x 50 mm
slope error tangential: 0.25 $\mu\text{rad rms}$
slope error sagittal: 5 $\mu\text{rad rms}$
HSF Roughness < 2 \AA rms
stripe 1: Ru / C
Stripe 2: W / B₄C
Period: 40 / 26 \AA
Density: Ru~10.5 g/cm³ C~2.2 g/cm³
W~17.5 g/cm³ B₄C~2.2 g/cm³
Interface Roughness 3 \AA rms



Synchrotron mirrors for PAL in Korea 300 mm in length (Operating energy range 10 – 80 KeV)

Montel multilayer mirror for tests and characterization and for the analyzer system of the ultra-high-resolution IXS spectrometer

2D parallel beam multilayer optic for 9100 eV
Montel part: 120 x 20 x 20 mm ± 0.1 mm (L x W x H)
Acceptance angle: > 10 mrad x 10 mrad
Meridional slope error: ≤ 2 arcsec (10 μrad) rms
Microroughness: ≤ 0.2 nm (rms) HSF
Coating: W/C Multilayer




Diamond Light Source Ltd
Beamline: B16 at DLS, 34-ID at APS
Contact: John Shutter, John.Sutter@Diamond.ac.uk


Tests and characterization of a laterally graded multilayer Montel mirror
J. Synchrotron Rad. (2014). 21, 16–23

Brookhaven National Laboratory, NSLS II
Beamline: 10-ID Inelastic X-ray Scattering
Contact: Yong Cai, cai@bnl.gov

Synchrotron X-ray tests of an L-shaped laterally graded multilayer mirror for the analyzer system of the ultra-high-resolution IXS spectrometer at NSLS-II
J. Synchrotron Rad. (2011). 18, 862–870




J. Synchrotron Rad. (2014). 21, 16–23




J. Synchrotron Rad. (2011). 18, 862–870

Collimating Montel mirror as part of a multi-crystal analyzer system for resonant inelastic X-ray scattering

2D parallel beam multilayer optic for 11215 eV
Montel part: 150 x 7 x 10 mm ± 0.1 mm (L x W x H)
Acceptance angle: > 14 mrad x 14 mrad
Meridional slope error: ≤ 2 arcsec (10 μrad) rms
Microroughness: ≤ 2.0 nm (rms) HSF
Coating: Ru/C Multilayer



Argonne National Lab, Advanced Photon Source
Beamline: 27-ID
Contact: Jung Ho Kim, jhkim@aps.anl.gov




J. Synchrotron Rad. (2016). 23, 880–886

Performance of a collimating L-shaped laterally graded multilayer mirror for the IXS analyzer system at NSLS-II

2D parallel beam multilayer optic for 9130 eV
Montel part: 100 x 4 x 6 mm ± 0.1 mm (L x W x H)
Acceptance angle: > 10 mrad x 10 mrad
Meridional slope error: ≤ 2 arcsec (10 μrad) rms
Microroughness: ≤ 0.2 nm (rms) HSF
Coating: W/C Multilayer




Brookhaven National Laboratory, NSLS II
Beamline: 30-ID beamline at APS
Contact: Yong Cai, cai@bnl.gov



J. Synchrotron Rad. (2014). 21, 473–478

Comparison of Scatterless Pinholes SCATEX and Scatterless Slits 2.0

The measurements were performed by C. Gollwitzer at the PTB four-crystal monochromator beamline at BESSY II at 8 keV with typical photon fluxes of ~10¹⁰ ph/s.




Images of the parasitic aperture scattering at 8 keV with the test apertures being the beam defining element. No scatter guard inserted. The downstream photon flux was the same (variation < 1%) for all compared test apertures.

Deduced scattering intensity vs. q -plots (360°-integration) for the various tested apertures.

- up to 4 times less parasitic aperture scattering compared to Scatterless Slits 2.0
- faster aperture scattering decay below the background at considerably smaller q -values

Comparison of SCATEX-Ta Pinholes and Tungsten Slits

The measurements were performed at 13 keV at the Nanofocus Endstation P03 beamline at PETRA III with typical photon fluxes of 10¹¹-10¹² ph/s by C. Krywka.



Detector images of the parasitic aperture scattering at 13 keV. In the standard beamline setup S5 denotes the position of the beam defining aperture and S6 the position of the antiscatter aperture.

Scattering intensity vs. q -plot. The data is normalized to the number of summed up pixel. Various apertures were tested at position S5 (beam definition) and S6 (scatter guard).

- a single SCATEX-Ta pinhole replaces both beam defining slit S5 and antiscatter slit S6
- the beam-defining SCATEX-Ta aperture can be positioned closer to the sample
- one order of magnitude less parasitic aperture scattering with SCATEX pinholes