

Trends on Multilayer X-ray Optics and Scatterless Apertures for X-ray Analytical Equipment

J. Graf¹, A. Kleine¹, S. Freisz², V. Smith², A. Beerlink¹, J. Wiesmann¹, C. Michaelsen¹

¹ Incoatec GmbH, Geesthacht, Germany, ² Bruker AXS GmbH, Karlsruhe, Germany



Introduction

X-ray optical components, such as multilayer mirrors or scatterless apertures, are used as beam conditioning devices in nearly all state-of-the-art X-ray analytical equipment, either in the home lab or at synchrotron beamlines. In this contribution, we will give an overview of current developments in multilayer optics and scatterless beam components, and show their benefit in combination with high-brightness microfocus X-ray sources for typical applications in small angle scattering and crystallography.

Multilayer Mirrors for Home Lab sources

Montel Multilayer Optics for 2D Beam Shaping

Montel optics are two mirrors mounted side-by-side in an L-shape enabling a 2-dimensional beam shaping. The mirrors comprise multilayer coatings that are typically deposited with a precision within $\pm 1\%$ of the d -spacing by physical vapor deposition techniques. For high-brightness laboratory sources, such as the latest Incoatec Microfocus Source μ S or the recently introduced liquid metal jet X-ray source, optics with similarly low shape errors are used delivering a well-shaped beam with a Gaussian-like intensity profile. A Montel optics with two elliptically shaped mirrors is point focusing (for e.g. single crystal diffraction), whereas two parabolic mirrors enable a collimated beam (for e.g. small angle scattering). Nowadays, Montel optics are also used at synchrotrons, where they substitute the KB (Kirkpatrick-Baez) mirrors achieving a more compact design.



Typical layer pair thickness accuracy for graded multilayers determined by X-ray reflectometry: we achieve a d spacing gradient accuracy of $\pm 1\%$; in the case of a uniform d -spacing we reach an accuracy $< \pm 0.1\%$.

Montel optics (top) and corresponding optical scheme (bottom) of a focusing multilayer mirror.


Multilayer Optics and Modern High-Brightness X-ray Sources

Generator	Incoatec μ S ^{High Brilliance}	Bruker TXS HB	Excillum METALJET
Optics	HELIOS MX	HELIOS MX	HELIOS MX GA
Power	50 W, 1-P	2500 W, 3-P	200 W, 1-P
Cooling	Air	Water / Air	Air
Wavelength [Å]	1.54	1.54	1.35
Spectral purity [% K α]	99.9	99.9	99.9
Beam diameter FWHM [μ m]	95	160	70
Divergence [mrad]	7.6	7.6	7.6
Intensity [ph/s/mm ²]	$> 2 \times 10^{10}$	$> 1 \times 10^{11}$	$> 4 \times 10^{11}$


Performance comparison of high-brightness microfocus X-ray sources and focusing multilayer mirrors for small molecule and protein crystallography.

Protein Crystallography and Screening

Sample: HIV Protease Complex	
$a = 46.38 \text{ \AA}$, $b = 57.86 \text{ \AA}$, $c = 84.87 \text{ \AA}$; P2 ₁ 2 ₁ 2 ₁ ;	
$T = 100 \text{ K}$; 99 amino acids	
crystal size	0.16 x 0.10 x 0.02 mm ³
exposure time	30 s/0.5°
total time	2 h
resolution	34 – 1.90 Å (2.00 – 1.90)
$\langle 1/\sigma \rangle$	14.5 (3.6)
$\langle \text{multiplicity} \rangle$	5.4 (5.4)
$\langle \text{completeness} \rangle$	99.9 % (100 %)
R_{int}	0.0855 (0.4331)



Data statistics and typical diffraction pattern of a small protein crystal, measured with a Bruker AXS D8 VENTURE equipped with the Excillum METALJET X-ray source.




Bruker D8 VENTURE equipped with an μ S^{High Brilliance} and an automated ISX stage for in-situ crystal screening and data collection in crystallisation plates at RT; typical diffraction pattern of a Thaumatin crystal.

SCATEX Pinholes for Home Lab Sources and Synchrotrons

Comparison of Tungsten Slits and SCATEX-Ta Pinholes

The measurements were performed at 13 keV at the Nanofocus Endstation P03 beamline at PETRA III with typical photon fluxes of 10^{11} - 10^{12} ph/s by C. Krywka.




Images of the parasitic aperture scattering at 13 keV. In the standard beamline setup S5 denotes the position of the beam defining aperture and S6 the position of the antiscatter aperture.

Scattering intensity vs. q -plot. The data is normalized to the number of summed up pixel. Various apertures were tested at position S5 (beam defining) and S6 (scatter guard).

- a single SCATEX-Ta pinhole replaces both beam defining slit S5 and antiscatter slit S6
- the beam-defining SCATEX-Ta aperture can be positioned closer to the sample
- one order of magnitude less parasitic aperture scattering with SCATEX pinholes
- pinhole sizes down to 10 - 20 μ m possible

Comparison of Scatterless Slits 2.0 and SCATEX Pinholes

The measurements were performed by C. Gollwitzer at the PTB four-crystal monochromator beamline at BESSY II at 8 keV with typical photon fluxes of $\sim 10^{10}$ ph/s.




Images of the parasitic aperture scattering at 8 keV with the test apertures being the beam defining element. No scatter guard inserted. The downstream photon flux was the same (variation $< 1\%$) for all compared test apertures.

Deduced scattering intensity vs. q -plots (360°-integration) for the various tested apertures.

New generation SCATEX pinholes

- up to 4 times less parasitic aperture scattering compared to Scatterless Slits 2.0
- up to 19 times less parasitic aperture scattering compared to old generation SCATEX Pinholes
- faster aperture scattering decay below the background at considerably smaller q -values


Comparison of a SCATEX 2-Pinhole Setup and a Standard 3-Pinhole Setup



SAXS Image of a thin fiber of a rat tail tendon, measured on a Bruker NANOSTAR with an μ S.

Advantages of a SCATEX 2-Pinhole Setup


- higher flux and smaller q_{min} possible due to a larger beam defining pinhole and a smaller beamstop
- faster data acquisition possible
- smaller footprint due to less pinholes and shorter beam path



Scattering intensity vs. q -plot, measured with a 3-pinhole high resolution NANOSTAR and a modified 2-pinhole NANOSTAR equipped with SCATEX pinholes. With a similar resolution the SCATEX setup gives a considerably higher scattering intensity.

Upgrading Existing Diffractometers

Incoatec supports the adaptation and full integration of our mirrors, SCATEX pinholes and μ S sources on almost all existing diffractometers.



μ S Upgrade on a mar-research 345 dtb (University of Basel, Switzerland)

μ S Adaptation on a HRXRD setup at synchrotron beamline (Petra III, DESY, Hamburg, Germany)

μ S and SCATEX Upgrade on a customized SAXS setup (University of Hamburg, Germany)